



Low-Temperature Thermal Properties of Carbon-Based Nanomaterials: Exploration of Graphene Oxide and Fullerene Composites



O. O. Romantsova^{1,2}, D. Szewczyk², A. Jeżowski², Y. V. Horbatenko¹,
A. I. Krivchikov^{1,2}, M. Vinnikov¹, S. Cherednichenko¹

¹B. Verkin Institute for Low Temperature Physics and Engineering of the NAS of Ukraine,
47 Nauky Ave., Kharkiv, 61103, Ukraine

²Institute of Low Temperature and Structure Research, PAS, 2 Okólna Str., 50-422, Wrocław, Poland

This study investigates the heat capacity of thermally reduced graphene oxide (TRGO-300) combined with fullerene C₆₀ (HyFn) within the temperature range of 2 K to 300 K, utilizing a Physical Property Measurement System (PPMS) with a Heat Capacity Module. In addition to these materials a reference tests with commercial samples from SigmaAldrich were performed both for pure fullerene C₆₀ and graphene oxide. The research aims to uncover how structural disorder influences thermal behavior, focusing on identifying anomalies such as the boson peak – a feature associated with vibrational density of states in disordered systems, crucial for understanding low-temperature thermal properties in nanostructured materials.

Samples:




Fig. 1 Examples of investigated materials on the sample platform (PPMS)




Fig. 2 Experimental data from TEM (20 nm).




Fig. 5. Experimental temperature dependence of heat capacity for TRGO with C₆₀, pristine TRGO, and pure C₆₀.




Fig. 8 Specific heat of the TRGO sample plotted as C_p/T^3 vs T

Raman spectra




Fig. 3 Raman-spectroscopy data for TRGO with C₆₀.

Experimental results:




Fig. 6. Specific heat of the TRGO with fullerene after deducting the contribution of tunnelling two-level systems plotted as C_p/T^3 vs T .




Fig. 9. Specific heat of the TRGO with C₆₀ sample plotted as C_p/T vs T^2 .

Measurement setup



Fig.4 The specific heat was determined using the method of thermal relaxation on a commercially available Physical Properties Measurement System (PPMS, Quantum Design Inc.)




Fig. 7. Specific heat of the C₆₀ sample plotted as C_p/T^3 vs T .




Fig. 10. Low-temperature specific heat data with fitting data using the soft potential model.

The main outcome of this study is the detection of a pronounced boson peak in TRGO–C₆₀ composites, revealed upon subtraction of the TLS (linear) contribution from the specific heat data. This study not only enhances understanding of heat transport mechanisms in carbon-based nanostructures but also provides valuable insights for applications in nanotechnology and cosmic science, where efficient thermal management is vital for performance and material stability. The findings underscore the importance of material synthesis methods in tailoring thermal properties, paving the way for future research on advanced carbon composites in extreme environments.

O. Romantsova gratefully acknowledges the support from the COST Action “Carbon molecular nanostructures in space (NanoSpace)” CA21126, supported by COST (European Cooperation in Science and Technology).